






LL Vernebler

CE 0123

| | | |
|-----------|----------------------------|-----|
| de | Gebrauchsanweisung | 1 |
| en | Instructions for use | 23 |
| da | Brugsanvisning..... | 45 |
| el | Οδηγίες χρήσεως | 65 |
| es | Instrucciones de uso..... | 87 |
| fi | Käyttöohje | 109 |
| fr | Notice d'utilisation..... | 129 |
| it | Istruzioni per l'uso | 151 |
| nl | Gebruiksaanwijzing | 173 |
| no | Bruksanvisning..... | 195 |
| pt | Manual de Instruções | 215 |
| sv | Bruksanvisning..... | 237 |
| tr | Kullanım kılavuzu | 257 |

A**B****C**

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Tel.: +49 (0)8151-279 279

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en Instructions for use

Read these instructions carefully and completely before using for the first time. Keep these instructions in a safe place.

Personal safety cannot be guaranteed if these instructions are not strictly followed.

Table of contents

| | | |
|----------|--|-----------|
| 1 | Components | 24 |
| 2 | Intended use | 24 |
| 3 | Information | 25 |
| 4 | Preparing to inhale | 26 |
| 4.1 | Assembling the nebuliser..... | 26 |
| 4.2 | Filling with medication..... | 27 |
| 5 | Inhalation | 28 |
| 5.1 | Important notes..... | 28 |
| 5.2 | Performing the inhalation..... | 28 |
| 6 | Hygienic re-use | 29 |
| 6.1 | General..... | 29 |
| 6.2 | Preparation..... | 30 |
| 6.3 | At home..... | 30 |
| 6.4 | Hospital and medical practice..... | 32 |
| 7 | Miscellaneous | 35 |
| 7.1 | Nozzle blocked..... | 35 |
| 7.2 | Material resistance..... | 35 |
| 7.3 | Storage and transportation conditions..... | 36 |
| 7.4 | Disposal..... | 36 |
| 8 | Technical data | 37 |
| 8.1 | General..... | 37 |
| 8.2 | Aerosol characteristics in accordance with DIN EN 13544-1 appendix CC..... | 37 |
| 9 | Replacement sets, replacement parts and accessories | 39 |

1 Components

Figure A:

Please refer to the following table to ensure that all illustrated components are included in the corresponding nebuliser kit package. If not, please notify the dealer from whom you obtained this PARI nebuliser immediately.

| | |
|-----------------------------|-----------------------------|
| (1) Inspiratory valve | (7) Nozzle cleaner |
| (2) Valve shim | (8) Nebuliser lower section |
| (3) Chamber | (9) Nozzle |
| (4) Nebuliser upper section | (10) Airflow control |
| (5) Medication cup | (11) Mouthpiece |
| (6) LL adapter | |

2 Intended use

The PARI LL nebuliser is a medication nebuliser which can be combined with the PARI compressor used or the PARI CENTRAL and the hospital's central compressed air supply and other PARI accessories as required to form a highly effective inhalation device for inhalation therapy of the lower respiratory tract. It is only suitable for temporary¹⁾ oral inhalation of liquid medicine prescribed or recommended by a doctor. The PARI LL nebuliser is capable of nebulising all commercially available medications that are currently approved for inhalation therapy.

The PARI LL nebuliser can be operated with **all** PARI compressors that are available on the market (for example the PARI BOY S, PARI BOY SX, PARI BOY N, PARI MASTER etc.) or with a clinic's central compressed air supply (see Chapter 8.1 "General"). **Please read and follow the instructions for use with your PARI compressor and/or PARI accessory.** PARI GmbH cannot be held liable for damage that is caused by using PARI nebulisers with other compressors.

1) Each application lasts less than 1 hour, in this case max. 20 minutes.

The manufacturer will be able to tell you whether the PARI LL nebuliser can be used with new PARI compressors that are introduced in the future. This information will also appear in the instructions for use of the compressor.

3 Information

The PARI LL nebuliser is a robust and durable nebuliser with an integrated interrupter button. This enables intermittent nebulising, which ensures optimum administration of the medicine.

The nebuliser can be combined with numerous practical PARI accessories and can thus be tailored in a variety of ways to the individual requirements of the patient.

- Only use medication prescribed or recommended by your doctor.
- You must conscientiously follow the general hygiene steps for using your nebuliser, e.g., thorough washing of hands or disinfection if required.
- The PARI LL nebuliser is only suitable for patients who are able to breathe spontaneously.
- In patients who have undergone tracheostomy, the PARI TRACHEO SET has to be used.

- Please follow the hygiene regulations before the device is used for the first time, before it is re-used, and after any long break in use.
- Ensure that the nebuliser was cleaned, disinfected, and, if appropriate, sterilised after the last treatment (see Chapter 6 "Hygienic re-use").

4.1 Assembling the nebuliser

Figure A:

The individual nebuliser parts are assembled as follows:

- Push the airflow control (10) as far as it will go onto the nozzle (9) and then push both parts into the medication cup from below (5).
- Insert the medication cup in the nebuliser lower section (8).
- Screw the nebuliser upper section (4) onto the nebuliser lower section.
- Slide the chamber (3) into the nebuliser upper section.
- Put the inspiratory valve (1) with the valve shim (2) on the chamber.
- Attach the LL adapter (6) to the nebuliser lower section (8) and connect the connection tubing (12) to the LL adapter and to the compressor or the PARI CENTRAL.
- Attach the mouthpiece (11) or the child or adult mask (13).

Safety instruction for PARI Child and Adult masks soft:

It is essential to ensure that the expiratory valve plate on your mask is pressed out before each use (see Figure C).

Inhalation with the mouthpiece is the most effective form of inhalation, since it results in minimum medication loss on the way to the lung. If it is not possible to use a mouthpiece, or if the patient prefers to inhale through a mask, PARI can provide the right accessories.

4.2 Filling with medication

Figure B:

- Remove the inspiratory valve (1) from the nebuliser.
- Pour the quantity of medication prescribed by your doctor into the nebuliser chamber (3).

Ensure that the medication is filled no higher than the upper scale marking (max. level 8ml). If you overfill the nebuliser and the medication leaks out of the bottom, please empty and clean (see Chapter 6 "Hygienic re-use").

Then fill the nebuliser with the medication again.

- Put the inspiratory valve (1) with the valve shim (2) on the chamber.

5.1 Important notes

- Before the patient starts inhaling, ensure that all parts are firmly connected to each other.
 - Ensure that the airflow control (10) is fully seated on the nozzle (9).
 - Ensure that the air outlets are connected firmly to the air connection on the PARI compressor or the PARI CENTRAL and on the nebuliser.
- Nebulisation can be adversely affected by an incorrectly assembled nebuliser.

Babies, infants, children, and anyone who requires assistance must be supervised constantly by an adult during inhalation therapy. Otherwise, safe and effective treatment cannot be guaranteed.

Caution: Contains small parts which can be swallowed!

5.2 Performing the inhalation

- Sit in an upright position and relax (see Figure D).
- Switch the compressor on or attach the connection tubing to the PARI CENTRAL and insert it in the compressed air outlet in the hospital wall provided for this purpose.
- Take the mouthpiece between your teeth and surround it with your lips (see Figure E) or hold the mask firmly over the mouth and nose.
- Press the interrupter button (see Figure F).
- Breathe in through the mouthpiece or mask as slowly and deeply as possible (see Figure G).
- Release the interrupter button and breathe in a relaxed manner through the mouthpiece or mask. The exhaled air should pass through the expiratory valve in the mouthpiece or the valve plate in the mask (expiratory valve will open, see Figure H).
- Continue breathing in and out as described until a change in the sound of the nebuliser signals that the medicine has been used up.
- Check whether a clearly visible aerosol is still flowing from the nebuliser (press interrupter button).
- Stop treatment as soon as the aerosol emerges irregularly.
- For continuous nebulising, press the interrupter switch and turn it gently clockwise to lock it in place.

6 Hygienic re-use

6.1 General

In order to prevent a health risk, e.g., infection caused by a contaminated nebuliser, it is imperative that the following hygiene regulations are followed.

The PARI LL nebuliser is designed for multiple uses. Please be aware that different requirements regarding hygienic re-use are in effect for different fields of application:

- At home (single patient use):
The nebuliser must be cleaned and disinfected after each use.
- Hospital / medical practice (possible change in patient):
The nebuliser must be cleaned, disinfected and sterilised after each use.
- Additional requirements concerning the necessary hygienic preparations (hand care, handling of medications or inhalation solutions) in high-risk groups (e.g., patients with cystic fibrosis) can be obtained from the respective self-help groups.
- Proof of suitability of the nebulisers for effective cleaning, disinfection and sterilisation has been obtained through an independent testing laboratory using the recommended procedures stated in Chapter 6.4.1 "Cleaning/disinfection" (test documentation available for inspection on request). Using the named alternatives remains the responsibility of the user.
- Please ensure adequate drying after each cleaning, disinfection and/or sterilisation. Condensation or residual wetness can present an increased risk through the growth of bacteria.
- Please check the parts of your nebuliser regularly and replace any defective (broken, misshapen, discoloured) parts.

The PARI Child and Adult soft masks cannot be sterilised.

Replace the nebuliser at least every two years. PARI GmbH recommends that certain components of the nebuliser (nozzle, connection tubing etc.) be replaced every year. For this purpose, PARI GmbH can offer you the useful PARI Year Pack (see Chapter 9 "Replacement sets, replacement parts and accessories"). **Please read and follow the hygiene instructions in the instructions for use of the PARI compressor.**

6.2 Preparation

Directly after each treatment, all parts of the PARI LL nebuliser must be cleaned of medication residue and contamination. To do this, dismantle the nebuliser into its individual parts (as described in the following).

Figure A:

- Detach the mouthpiece (11) or the mask (13), the inspiratory valve (1) and the chamber (3) from the nebuliser.
- Remove the valve shim (2) from the inspiratory valve.
- Unscrew the nebuliser upper section (4) from the nebuliser lower section (8).
- Push the medication cup (5) up as far as the windows and then pull it out of the nebuliser lower section.
- Press the nozzle (9) together with the airflow control (10) down and out of the medication cup.
- Pull the airflow control up and off from the nozzle.
- Turn the interrupter button clockwise as far as it will go.

For hygienic re-use at home please continue reading the following section.

For hygienic re-use in the hospital or medical practice, please read Chapter 6.4 "Hospital and medical practice".

6.3 At home

6.3.1 Cleaning

- Remove the tubing from the nebuliser and dismantle the nebuliser (see Chapter 6.2 "Preparation").

Refer to Section 6.3.5 for care of the connection tubing.

Recommended method: Use warm tap water

- Clean all parts of the nebuliser thoroughly for 5 minutes with warm tap water (approx. 40°C) and a little dishwashing liquid (dosage according to the manufacturer's recommendations).
- Then rinse all parts thoroughly under warm running water (approx. 40°C without dishwashing liquid).
- Drain off well.

You can remove excess water more quickly by shaking the parts.

Possible alternative: Wash in the dishwasher

- Clean the nebuliser in your dishwasher (not together with dirty crockery).
- Distribute the nebuliser parts in the cutlery basket.
- Select a programme of at least 50°C.

6.3.2 Disinfection

Disinfect the dismantled nebuliser after cleaning.

Recommended method: Using a standard disinfectant for baby bottles

Disinfection must last at least 15 minutes. As regards disinfection and the quantity of water necessary make sure you observe the instructions for use of the relevant disinfectant. Always check carefully that it is clean and fully functional.

Possible alternative: In boiling water

- Place the individual parts of your nebuliser in boiling water for at least 15 minutes.
Use a clean saucepan. If your local tap water is hard, please use low-lime or distilled water.

Ensure that the water level in the pot remains high enough to prevent the nebuliser parts from coming into contact with the hot base of the pot.

6.3.3 Drying

- Place the nebuliser parts on a dry, clean and absorbent surface and let them dry completely (at least 4 hours).

6.3.4 Storage

- Between uses, and in particular during long breaks in treatment, wrap the nebuliser in a clean, lint-free cloth (e.g., a tea-towel) and keep it in a dry, dust-free environment.

6.3.5 Care of the connection tubing

- After each treatment, pull the tubing from the nebuliser and switch on your PARI compressor.
- Leave the compressor running until any condensation in the tubing has been removed by the airflow.
- Replace the connection tubing if it is contaminated (part no. 041B4591, see Chapter 9).

Please continue reading Chapter 7 "Miscellaneous".

6.4 Hospital and medical practice

When the nebuliser is used for several patients, it must be cleaned, disinfected and sterilised before each change in patient.

Please ensure that only adequately validated device-specific and product-specific procedures are used for cleaning, disinfection and sterilisation, and that the validated parameters are observed in each cycle.

The cleaning and disinfection procedures must be recognised as effective (e.g., in the list of tested and recognised disinfectants and procedures issued by the Robert Koch Institute/DGHM) and already validated in all cases.

If other procedures are used, their effectiveness must be confirmed during validation. CE-approved chemical disinfectants may also be used, provided they satisfy the requirements for material resistance (see Chapter 7.2 "Material resistance").

Please follow the instructions in Chapter 7.2 "Material resistance".

In addition, please follow the hygiene regulations of the hospital or medical practice.

6.4.1 Cleaning/disinfection

Clean and disinfect immediately after use. Ideally, a mechanical procedure (instrument dishwasher) should be used.

- Prepare for cleaning and disinfection as described in Chapter 6.2 "Preparation".

Any suitable cleaner/disinfectant can be used.

Recommended procedure: Thermal disinfection

- Place the nebuliser components in the instrument dishwasher.
- Select the 93°C programme (treatment time 10 minutes).

The effectiveness of this procedure has been confirmed using a G7836 CD disinfectant manufactured by Miele with neodisher®medizym cleaning agent manufactured by Chemische Fabrik Dr. Weigert, Hamburg.

or

Recommended procedure: Manual cleaning/disinfection

The effectiveness of this procedure has been confirmed using the Korsolex®-Endo-Cleaner cleaning agent in combination with the Korsolex®basic disinfectant manufactured by BODE CHEMIE HAMBURG.

- The components must be placed in a 0.5% solution of Korsolex®-Endo-Cleaner and approx. 50°C warm water, and cleaned thoroughly for 5 minutes.
- The components are then placed in a 4% disinfectant solution of Korsolex®basic for 15 minutes.
- Then thoroughly rinse the nebuliser components under hygienically tested water and leave them to dry completely on a clean, dry and absorbent surface (for at least 4 hours).

Possible alternative: Chemical-thermal disinfection

- Place the nebuliser components in the instrument dishwasher.
- Add a suitable cleaner/disinfectant and select the 60°C programme.

Storage

- Between uses, wrap the nebuliser in a clean, lint-free cloth and keep it in a dry, dust-free place protected from contamination.

6.4.2 Sterilisation

After cleaning/disinfection, pack the dismantled nebuliser in sterilisation packaging (disposable sterilisation packaging, e.g., foil/paper sterilisation bags). The sterilisation packaging must conform to DIN EN ISO 11607 and be suitable for steam sterilisation. Then sterilise the nebuliser using the following sterilisation procedure:

Recommended procedure: Steam sterilisation

Sterilisation temperature: 121°C (to be maintained for at least 20 min.) or 132°C/134°C (to be maintained for at least 3 min.) max. 137°C.

This procedure has been validated according to DIN EN ISO 17665 as well as DGKH guidelines for the validation and routine monitoring of sterilisation processes using damp heat for medical devices.

Storage

Keep the sterilised nebuliser in a dry, dust-free place protected from contamination.

6.4.3 Care of the connection tubing

- If condensation is visible in the tubing, pull the tubing off the nebuliser after treatment and switch on your PARI compressor.
- Leave the compressor running until any condensation in the tubing has been removed by the airflow.
- Replace the connection tubing in the case of a change of patients or if it becomes contaminated (part no. 041G4550, see Chapter 9).

7 Miscellaneous

7.1 Nozzle blocked

Should the nozzle be blocked, you can clear it yourself as follows.

Figure 1:

- Unscrew the nebuliser upper section (4) from the nebuliser lower section (8).
- Push the medication cup (5) up as far as the windows and then pull it out of the nebuliser lower section.
- Press the nozzle (9) together with the airflow control (10) down and out of the medication cup.
- Pull the airflow control up and off from the nozzle.
- Take the nozzle cleaner (7) out of the nebuliser lower section.
- Carefully poke through the three nozzle bore holes "A" to "C" (only from below).
- Re-assemble the nebuliser (see Chapter 4.1).

- Only use an original PARI nozzle cleaner!
- Use the nozzle cleaner very carefully, because the small nozzle bore holes can easily be damaged irreparably.
- If the nozzle becomes blocked repeatedly, replace it (part no. 019B2110).

7.2 Material resistance

Do not expose the nebuliser to temperatures above 137°C. In general, like all plastic parts, nebulisers are susceptible to wear. Nebulisers can be sterilised up to 300 times, but then they must be replaced. Even if the nebuliser has been sterilised fewer than 300 times, it must still be replaced at least every 2 years. PARI GmbH recommends that certain components of the nebuliser (nozzle, connection tubing etc.) be replaced every year. Note the following when selecting cleaners/disinfectants:

- In general, the aldehydic group of cleaners/disinfectants is suitable for cleaning and/or disinfecting the PARI LL nebuliser.
- The use of other cleaner or disinfectant groups has not been tested with regard to the material resistance of this nebuliser.

7.3 Storage and transportation conditions

Do not store in damp rooms (e.g., bathrooms) or transport together with damp objects.

Protect from continuous direct sunlight during storage and transportation.

In general, the nebuliser should be stored as described in Chapter 6

"Hygienic re-use", "Storage" between applications. See Chapter 6.3.4 (At home) or Chapter 6.4 (Hospital and medical practice).

7.4 Disposal

The PARI LL nebuliser may be disposed of with domestic waste unless the disposal regulations in the respective member countries provide otherwise (e.g., disposal by local authorities or dealers).

Materials recycling helps to reduce the consumption of raw materials and protect the environment.

8 Technical data

8.1 General


| | |
|-----------------------------------|---------------|
| Recommended nebuliser gases: | Air, oxygen |
| Minimum flow: | 3.5l/min |
| Corresponding operating pressure: | 0.8bar/80kPa |
| Maximum flow: | 6.0l/min. |
| Corresponding operating pressure: | 2.0bar/200kPa |
| Minimum fill: | 2ml |
| Maximum fill: | 8ml |

8.2 Aerosol characteristics in accordance with DIN EN 13544-1 appendix CC

8.2.1 Determination of particle size

Fill the nebuliser with 3ml of a 2.5% solution of sodium fluoride and then attach the nebuliser to a Marple cascade impactor and run for at least 3 minutes.

Cumulative mass distribution for the PARI LL nebuliser at the minimum flow (3.5l/min) and the maximum flow (6.0l/min):

(Averages from each of 2 separate measurement series with 3 nebuliser/ compressed air combinations)

8.2.2 Determination of the aerosol output and aerosol output rate

Aerosol output (aerosol quantity output by the nebuliser system for a specified fill volume)

Fill the nebuliser with 2ml of a 1% sodium fluoride solution. Attach the nebuliser to the PARI Compass Breath Simulator until the aerosol becomes irregular, then one more minute.

Aerosol output in [ml] for minimum flow (3.5l/min) and maximum flow (6.0l/min):

| | |
|----------------------|--------|
| Min. flow (3.5l/min) | 0.43ml |
| Max. flow (6.0l/min) | 0.50ml |

(Averages from each of 2 separate measurement series with 3 nebuliser/compressed air combinations)

Aerosol output rate (volume of aerosol that is output by the nebuliser system per unit of time)

Fill the nebuliser with 2ml of a 1% sodium fluoride solution. Attach the nebuliser to the PARI Compass Breath Simulator for at least 1 minute.





Aerosol output rate in [ml/min] for minimum flow (3.5l/min) and maximum flow (6.0l/min):





| | |
|----------------------|-------------|
| Min. flow (3.5l/min) | 0.113ml/min |
| Max. flow (6.0l/min) | 0.190ml/min |





(Averages from each of 2 separate measurement series with 3 nebuliser/compressed air combinations)





9 Replacement sets, replacement parts and accessories


For illustrations of replacement parts and accessories, refer to the fold-out pages in the brochure cover.

| | Description | Part no. |
|---|-------------------------|----------|
| Replacement sets: | | |
|  | PARI BOY Year Pack | 012G6300 |
| Replacement parts and accessories: | | |
|  | PARI LL nebuliser | 012G6202 |
|  | Nebuliser upper section | 012B1800 |
|  | Nebuliser lower section | 012B1200 |

| | Description | Part no. |
|--|---|-----------------|
|  | Medication cup | 012E1601 |
|  | Nozzle | 019B2110 |
|  | Inspiratory valve | 022B0510 |
|  | Universal mouthpiece (with expiratory valve) | 022E3050 |
|  | LL adapter | 012E1024 |
|  | Connection tubing for PARI compressors (1.2 m) | 041B4591 |

| | Description | Part no. |
|---|---|----------|
|  | Connection tubing (2.1 m; for use only with the hospital compressed air system) | 041G4550 |
|  | Nose clip | 041E3500 |
|  | Adult mask soft | 041G0740 |
| | Child mask soft "Spiggy" | 041G0741 |
|  | PARI SMARTMASK (inc. 45 degree knee) | 041G0730 |
|  | PARI Filter/Valve Set | 041G0500 |

| | Description | Part no. |
|--|---------------------------|----------|
|  | PARI PEP System I | 018G6100 |
|  | PARI PEP System II | 018G6200 |
|  | PARI PEP S System | 018G4000 |
|  | PARI THERM aerosol heater | 012G7110 |

| | Description | Part no. |
|---|------------------|----------|
|  | PARI TRACHEO SET | 025G0007 |

We reserve the right to make technical changes.

For all product information please contact our Service Center:

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