
**seca 959 /
seca 958**

D	<i>Bedienungsanleitung und Garantieerklärung</i>	3
GB	<i>Instruction manual and guarantee</i>	19
F	<i>Mode d'emploi et garantie</i>	35
I	<i>Manuale di istruzioni e garanzia</i>	51
E	<i>Manual de instrucciones y garantía</i>	67
DK	<i>Betjeningsvejledning og garantibevis</i>	83
S	<i>Bruksanvisning och garanti</i>	99
N	<i>Bruksanvisning og garantierklæring</i>	115
FIN	<i>Käyttöohje ja takuu</i>	131
NL	<i>Bedieningshandleiding en garantieverklaring</i>	147
P	<i>Instruções de utilização e declaração de garantia</i>	163
GR	<i>Οδηγίες χειρισμού και εγγύηση</i>	179

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Full certification

In seca products you are purchasing not only technology matured over a century, but also quality confirmed by authorities, the legal system and by institutes. seca products comply with European directives, standards and national laws. With seca, you are buying the future.

M

The products in these operating instructions comply with German legislation governing medical products and therefore with Directive 93/42/EEC of the Council of the European Community which is embodied in national laws throughout Europe.

Scales carrying this symbol are approved for medical purposes in the European Community. seca scales with this symbol meet the stringent quality and technical requirements placed on scales which can be calibrated.

Scales carrying this symbol are calibrated in accordance with Precision Class III of EC Directive 90/384/EEC.

Products carrying this symbol comply with the following directives and standards:

1. Directive 90/384/EEC governing non-automatic weighing instruments
2. Directive 93/42/EEC governing medical devices
3. DIN EN 45501 governing metrological aspects of non-automatic scales

seca's professionalism has also been officially recognized. The certificate from TÜV Product Service, the body responsible for medical devices, confirms that as a manufacturer of medical devices, seca rigorously complies with stringent legal requirements. seca's quality assurance system covers the development, production, sales and service of medical scales and measuring systems and ergometry.

seca helps the environment.

We are concerned to save natural resources. We therefore strive to save packaging material where meaningful. What is left over can conveniently be disposed of locally via the German Dual System (recycling programme).

1. Congratulations!

By purchasing the **seca 958/959** electronic chair scale, you have acquired a highly-accurate and sturdy piece of equipment.

seca has been putting its experience at the service of health for over 150 years now, and as market leader in many countries of the world, is always setting new standards with its innovative developments for weighing and measuring.

In accordance with national regulations, the **seca 958/959** wheelchair scale with standing aid is used primarily in hospitals, doctors' surgeries and in-patient nursing facilities. The scale is a class I medical product and is approved to CE class III (seca 959)/III (seca 958).

Weight is determined within a few seconds.

The **seca 958/959** can be moved on castors and its low power consumption combined with its rechargeable batteries mean it will provide long, mobile service.


seca 958/959 chair scales are made for people who are not especially steady on their feet. To give you a safe feeling while weighing, the chair scale has a solid tubular frame and a firmly bolted seat. The free-running casters make the scale very easy to move. The wheel-locking lever, together with the swiveling armrests and foot supports, makes it easier to sit in safety.

The scale is of a very sturdy construction and will give you long and loyal service. It is easy to operate and the large display is easy to read.

2. Safety

Before using the new scale, please take a little time to read the following safety instructions.

Warnings

- When adjusting the height of the screw feet of the chair scale, ensure that the screw is screwed at least 15 mm into the frame. If you screw the screw foot less than 15 mm into the frame, there is a risk of the screw foot breaking off and

the occupant of the chair scale being injured.

- Never use the **seca 958/959** chair scale for transport purposes. Do not use **seca 958/959** chair scale to transport people or objects. If you transport people or objects using the chair scale, there is an acute risk of injury.
- Never place the **seca 958/959** chair scale on sloping surfaces. The brakes on the casters could be released.
- Prevent the scale from rolling away accidentally. Use the wheel brake when the patient is sitting down or getting up and assist people who are unsteady on their feet.
- When using the scale with a mains unit, ensure that the supply cable is routed in such a way as to exclude any type of tripping hazard.

Notes

- Follow the safety instructions in the instruction manual.
 - Keep the operating manual and the declaration of conformity contained therein in a safe place.
 - Have scale serviced and re-calibrated on a regular basis (see “Servicing/recalibration” on page 31).
 - Be sure to use only the seca mains unit
- supplied. Before plugging the mains unit into a socket, check that the mains voltage indicated on the mains unit matches local mains voltage.
- Ensure that the cables are not crushed and do not come into contact with hot objects.
 - Use only the type of battery stated (see “Spare parts and accessories” on page 32).

3. Before you start...

Unpacking


The scope of supply includes:

- Chair scale comes complete or partly assembled
- 1 seca mains unit
- 1 rechargeable battery pack
- 1 set of operating instructions


Assembling

This chair scale is available either fully or partly assembled. To assemble the partly-assembled model, please proceed as follows:


- Locate the handles on the journals provided and bolt in place using the two M6 x 16 bolts, following the diagrams.
- Thread the mounting bolts for the display unit through the bores provided in the handles and screw the display tightly in place using two M4 nuts.


- Attach the weighing cell cable using the enclosed plastic cable clip by pushing this into the bores on the handles.

- Fit the seat to the frame using the four self-tapping plastic screws. Check that the seat is screwed tightly to the frame with no play at all.

- Push the hinge bolts for the armrests through the hinge bushes on the seat frame. Put on the washer and cap nut (M10) and screw them tight.


Power supply

The scale is powered by a mains unit and by rechargeable batteries.

Connecting the rechargeable battery pack:

- Open the battery compartment in the base of the display housing.
- Connect the plug contacts to the rechargeable battery pack.
- Close the lid of the battery compartment again.


The LC display and the evaluation electronics have a low power consumption – 1,000 weighing operations can be performed with a fully charged battery.

Once battery voltage is no longer adequate, the battery symbol  first appears in the display to indicate that battery power is dropping.

If the batteries are not charged, battery voltage will ultimately not be adequate for error-free measuring. In this case, **bAtE** appears in the display and weighing is no longer possible.



The rechargeable batteries are charged by the mains unit supplied. For this purpose, the mains unit needs to be connected to the mains for at least 24 hours.

Please use only the mains unit supplied or mains units from the seca range of accessories. Before plugging the unit into the socket, check whether the mains voltage marked on the mains unit matches local mains voltage.

4. Setting up and aligning the scale

- Place the scale on a level surface.
- Align the scale using the two bottom knurled nuts (1 and 2) and the two adjustable feet (3 and 4) so that the display in the spirit level (5) is central. This setting can be secured using the top lock-nuts (6, 7, 8 and 9).


spirit level

5. Operation

Controls and displays

Controls on the display

Switch scale on and off

press briefly → set display to zero, activate/deactivate tare function
in function menu → increase value

press briefly → activate/deactivate hold function
in function menu → decrease value


Activating the BMI function, Confirmation


Other controls

The armrests on the chair scale can be swivelled backwards independently. This is useful because it makes it easier for the person being weighed to sit down.


- Ensure that you do not trap any part of the body or pieces of clothing when folding the rests down.

- The rear casters can be locked if you push down the two locking levers on the casters. This device should be used when the patient is sitting down on the scale or getting up from it.

- Once the person who is to be weighed is seated, the foot-rests are pulled out and tilted downwards. The patient can now place his or her feet on the foot-rests.

Checking that the contents of the calibration counter are correct

This seca scale is calibrated. Calibration may be performed only by authorised centres. To ensure this, the scale is fitted with a calibration counter which records any change to data relevant to calibration. If you want to check whether the scale is properly calibrated, proceed as follows:

- If necessary, switch off the scale.
- Keep any key depressed and start the scale.

The current contents of the calibration counter are shown on the display for a few seconds.

- Compare the output contents of the calibration counter with the number shown on the calibration counter sticker.

For valid calibration, both numbers must match. If sticker and calibration counter do not match, then re-calibration is required. Contact your service agent or seca Customer Service.

START

2

seca
2

Both numbers must match.

If a supplementary calibration is necessary, then the supplementary calibration mark illustrated opposite is used to record the calibration counter position instead of the seca calibration mark shown above. This mark is secured with an additional seal by the person authorised to perform the supplementary calibration. The supplementary calibration mark is available from seca Customer Service under the number 14-05-01-886.

Correct weighing

- Press the Start key with no load on the scale.
SECA, **88888** and **00** appear consecutively in the display. The scale is then automatically set to zero and ready to use.
- Get the patient to sit on the scale.
- The patient's weight can only be determined once he or she is no longer in contact with the floor. Be sure, therefore, that the patient's feet are on the foot rest throughout the weighing operation.
- Read the weighing result off the digital display.
- If a load greater than 200 kg is placed on the scale, **STOP** appears in the display.
- To switch off the scale, press the Start key again.

Automatic switch-off/economy mode

The scale switches off automatically after 5 minutes.

Taring weight display (TARE)

The Tare function is very useful if an additional weight (a blanket, for example) is to be ignored for weighing purposes.

Proceed as follows:

- Press the Start key with no load on the scale.

- First place the additional weight (a blanket, for example) on the scale and activate the Tare function by pressing the **TARE** key.

The display is reset to zero. **0.0** appears in the display. The note "NET" appears.

- Now put the patient and the additional weight on the scale.

The scales calculate the patient's weight without the additional weight. You can now perform as many weighing operations as you like.

- Press the key again to deactivate the Tare function; the scale is then back in normal weighing mode and the tare value has been erased.
- The tare value is likewise erased by the scale being switched off.

Store weight value (HOLD)

You can also continue displaying the weight determined once the load has been removed. This makes it possible to look after a patient first, before having to note the weight:

- Press the **HOLD** key whilst the load is on the scale.

"HOLD" appears in the digital display. The value is stored until the scale switches off automatically.

The triangle with the exclamation mark indicates that this frozen display is an additional facility which is not subject to calibration.

- The HOLD function can be switched off by pressing the **HOLD** key again.

0.0 appears in the display if the scale has no load on it, otherwise the current weight. The "HOLD" display vanishes.

Note:

The HOLD function can also be used if the TARE function has been executed beforehand.

Body Mass Index (BMI)

Body Mass Index puts height and weight in relationship to one another, providing more accurate information than ideal Broca weight, for example. A tolerance range is given which is considered ideal for health.

- Press the Start key with no load on the scale.

SECR, **88888** and **00** appear consecutively in the display. The scale is then automatically set to zero and ready for operation.

- Press the **FUNC** key.

On the display, you will see the last height entered.

- You can use the arrow keys to adjust the value in increments according to the graduations displayed.

- Once the correct value is set, press the **FUNC** key again. The BMI function is now activated.

- Get the patient to sit on the scale.
- Read off the BMI on the digital display and compare it with the categories on page 29.

- Press an arrow key if you would like to see the weight briefly.

- To return to normal weighing mode, press the **FUNC** key.

- To switch off the scale, press the Start key again.

Evaluating Body Mass Index for adults:

Compare the value determined with the categories quoted below which correspond to those used by WHO, 2000.

- **BMI below 18.5:**

Warning, this patient weighs too little. There could be a tendency to anorexia. An increase in weight is recommended to improve well-being and performance. If in doubt, consult a specialist.

- **BMI between 18.5 and 24.9:**

The patient is of a normal weight. He can stay as he is.

- **BMI between 25 and 30 (pre-obese):**

The patient is slightly to moderately overweight. He should cut his weight if he is already ill (e.g. diabetes, hypertension, gout, fat metabolism disorders).

- **BMI over 30:**

Weight reduction is essential. Metabolism, cardiovascular system and bones are all under strain. A consistent diet, plenty of movement and behavioural training are all recommended. If in doubt, consult a specialist.

6. Cleaning

Clean the seat and the housing as required using a domestic cleaning agent or commercially-available disinfectant. Follow the manufacturer's instructions.

Under no circumstances use abrasive or acid cleaners, white spirit, benzene or the like for cleaning. Such substances can damage the high-quality surfaces.

7. What do I do if...

... no weight display comes on?

- Is the scale switched on?
- Check the mains connection and the connection of the battery pack.
- Switch the scale off for a few seconds. This can be done by opening the plug connected to the rechargeable battery block and unplugging the mains supply from the scale.

... 0.0 does not appear before weighing?

- Switch off the scale with the Start key and start the scale again – there must not be any load on the scale – and only its feet should be in contact with the floor.

... one segment is illuminated constantly or not at all?

- The relevant segment has a fault. Inform Servicing.

... the display shows ?

- Battery voltage is dropping. Charge the rechargeable batteries in the next few days.

... bAtt appears in the display?

- The rechargeable battery is so severely discharged that weighing is no longer possible. Recharge the battery.

... 5tOP appears in the display?

- Maximum load has been exceeded.

... the display flashes?

- If you have just activated the Hold or Tare function, please wait until the scale has detected a steady weight value. The display will then stop flashing.
- If there has been no activation of a function beforehand, remove the load from the scale and wait until 0.0 is displayed, then weigh again.

... the display tEMP appears?

- The ambient temperature of the scale is too high or too low. Place the scale in an ambient temperature between +10 °C and +40 °C. Wait about 15 minutes for the scale to adapt to the ambient temperature and then weigh again.

... the display E and a number appear?

- Switch off the scale with the Start key and start the scale again. The scale will then work normally again.

If this is not the case, disconnect the power supply by briefly removing the batteries. If this measure is equally unsuccessful, inform the Service department.

8. Servicing/recalibration

Have the scale recalibrated by authorised persons in accordance with national legal requirements. The year of initial calibration can be found behind the CE symbol on the rating plate above the number of the appointed office, 0109 (Hessen Weights & Measures Office). Recalibration is essential if one or more security marks are damaged or the calibration counter does not match the number on the applicable calibration counter mark (cf. page 26).

We recommend having your scale serviced by your local service agent before having it recalibrated. seca Customer Service will be pleased to give you further assistance.

9. Technical data

Maximum load	200 kg
Minimum load	1 kg
Graduations	
seca 959:	
0 up to 150 kg:	50 g
150 up to 200 kg:	100 g
seca 958:	
0 up to 100 kg:	100 g
100 up to 200 kg:	200 g
Dimensions	
height:	920 mm
width:	565 mm
depth:	810 mm
Weight	approx. 24.5 kg
Temperature range	+10 °C to +40 °C
Power supply	rechargeable batteries/mains unit
Power consumption	approx. 25 mA
Calibration as per Directive 90/384/EEC	
seca 959:	approved for medical use, Class III
seca 958:	approved for medical use, Class III
Medical product as per Directive 93/42/EEC	
EN 60 601-1:	Class I
Insulated equipment, Protection Class II:	
medical electrical equipment type B:	

Accuracy at first calibration

seca 959:

to 25 kg:	± 25 g
25 kg to 100 kg:	± 50 g
100 kg to 150 kg:	± 75 g
150 kg to 200 kg:	± 100 g

seca 958:

to 5 kg:	± 50 g
5 kg to 20 kg:	± 100 g
20 kg to 100 kg:	± 150 g
100 kg to 200 kg:	± 300 g

Functions

Tare, Hold, BMI

10. Spare parts and accessories

Power supply

Model dependent mains unit

Plug-in mains unit, Euro: 230V~ / 50Hz / 12V= / 130mA

seca no.: 68-32-10-252

Switchmode power adapter, 100-240V~ / 50-60Hz / 12V= / 0.5A

seca no. 68-32-10-265

Attention!

Use only seca original plug-in mains units with 9 V or controlled 12 V output voltage. Commercially-available mains units may supply a higher voltage than stated and thus lead to damage to the scale.

Model 471 power supply carrier

seca spare part no.: 471-00-00-009

Rechargeable batteries

seca spare part no.: 68-22-12-721

Recalibration

Supplementary calibration mark

seca spare part no.: 14-05-01-886

11. Disposal

Disposing of the device

Do not dispose of the device in domestic waste. The device must be disposed of properly as electronic waste. Follow the national regulations which apply in your case. For further information, contact our service department at: service@seca.com

Rechargeable batteries

Do not throw old rechargeable batteries into the domestic refuse! Ensure that rechargeable batteries are properly disposed of at local collection points.

12. Warranty

A two-year warranty from date of delivery applies to defects attributable to poor materials or workmanship. All moveable parts - batteries, cables, mains units, rechargeable batteries etc. - are excluded. Defects which come under warranty will be made good for the customer at no charge on production of the receipt. No further claims can be entertained. The costs of transport in both directions will be borne by the customer should the equipment be located anywhere other than the customer's premises. In the event of transport damage, claims under warranty can be

honoured only if the complete original packaging was used for any transport and the scale secured and attached in that packaging just as it was when originally packed. All the packaging should therefore be retained.

A claim under warranty will not be honoured if the equipment is opened by persons not expressly authorised by seca to do so.

We would ask our customers abroad to contact their local sales agent in the event of a warranty matter.

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